

IHE Integration Statement			
Vendor	Product Name	Version	Date
mediDOK Software Entwicklungsgesellschaft mbH	mediDOK_IHE_Bridge	2.1	4/15/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Consistent Time	Time Client	None
Mobile access to Health Documents	Document Consumer	None
Mobile access to Health Documents	Document Recipient	None
Mobile access to Health Documents	Document Recipient	XDS on FHIR
Mobile access to Health Documents	Document Responder	None
Mobile access to Health Documents	Document Responder	XDS on FHIR
Mobile access to Health Documents	Document Source	None
Cross-Enterprise Document Sharing	Document Consumer	Document Metadata Update
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Registry	Document Metadata Update
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V3)
Cross-Enterprise Document Sharing	Document Source	Document Replacement
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Integrated Document Source/Repository	Document Replacement

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
mediDOK Software Entwicklungsgesellschaft mbH	mediDOK_IHE_Bridge	2.1	4/15/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

<b>Integration Profile</b>	<b>Actor</b>	<b>Integration Profile Option</b>
Cross-Enterprise Document Sharing	Integrated Document Source/Repository	None
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Document Registry	None
Cross-Enterprise Document Sharing for Imaging	Document Repository	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	
<b>DICOM</b>	

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)